Externe Kontrollen *Serazym*[®] Stuhlteste

REF EK-017; EK-018; EK-020; EK-030; EK-039; EK-040;
EK-045; EK-061; EK-093; EK-106; EK-107; EK-114

Seramun Diagnostica GmbH · Spreenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
Telefon +49 (0) 33767 79110 · Fax +49 (0) 33767 79199 · info@seramun.com

Einführung

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Introduction

The External controls *Serazym*[®] Stool diagnostics are intended for use in external quality assessments for stool diagnostics of the *Serazym*[®] product line (Seramun Diagnostica GmbH) and for quality assurance of internal laboratory test processing.

The controls provide qualitative results (positive or negative).

Packaging Unit

Two External controls *Serazym*[®] Stool diagnostics are offered for every test parameter as a set with 400 µl stabilized parameter-specific antigen preparation each. The sets are shipped twice per year on predefined dates. The dates are announced by Virotech Diagnostics GmbH.

Preparation

Allow External controls *Serazym*[®] Stool diagnostics to reach room temperature prior to use and mix thoroughly.

External controls *Serazym*[®] Stool diagnostics are treated like sample material and are diluted with the respective sample diluent according to the manufacturer's instructions.

Note: The *Serazym*[®] Verotoxin 1+2 external controls do not contain viable bacteria. Therefore, they should not be pre-incubated in enrichment bouillon (potential false negative results by high dilution).

Intended Use

The External controls *Serazym*[®] Stool diagnostics can only be used in the following ELISA kits. There are no reference values for the use in ELISA kits of other manufactures.

Control-Set	Order-No.	ELISA Kit
Rotavirus-1 Rotavirus-2	EK-020	<i>Serazym</i> [®] Rotavirus
Adenovirus-1 Adenovirus-2	EK-017	<i>Serazym</i> [®] Adenovirus
Astrovirus-1 Astrovirus-2	EK-045	<i>Serazym</i> [®] Astrovirus
Norovirus-1 Norovirus-2	EK-061	<i>Serazym</i> [®] Norovirus
Campylobacter-1 Campylobacter-2	EK-093	<i>Serazym</i> [®] Campylobacter
Clostridium difficile GDH-1 Clostridium difficile GDH-2	EK-107	<i>Serazym</i> [®] Clostridium difficile GDH
Clostridium difficile Toxin A+B-1 Clostridium difficile Toxin A+B-2	EK-040	<i>Serazym</i> [®] Clostridium difficile Toxin A+B
Verotoxin 1+2-1 Verotoxin 1+2-2	EK-030	<i>Serazym</i> [®] Verotoxin 1+2
Giardia-1 Giardia-2	EK-106	<i>Serazym</i> [®] Giardia
Cryptosporidium parvum-1 Cryptosporidium parvum-2	EK-039	<i>Serazym</i> [®] Cryptosporidium parvum
Entamoeba histolytica-1 Entamoeba histolytica-2	EK-018	<i>Serazym</i> [®] Entamoeba histolytica
H. pylori 2nd Gen.-1 H. pylori 2nd Gen.-2	EK-114	<i>Serazym</i> [®] H. pylori 2nd Gen.

Assay Procedure

The controls are used like sample material without prior enrichment according to the manufacturer's instructions for the respective test kit.

Evaluation of Results

The evaluation of the External controls *Serazym*[®] Stool diagnostics is carried out according to the manufacturer's instructions of the respective test kit.

After the shelf life has expired, the qualitative results are published on the homepage of Virotech Diagnostics GmbH.


For this please enter the internet address <http://www.virotechdiagnostics.com/service.html> and choose External controls *Serazym*[®] Stool diagnostics. There you will find the certificate of analysis with the results.

Storage and Stability

The External controls *Serazym*[®] Stool diagnostics may only be used within the expiry date reported on its respective label (storage at 2...8°C).

History of Changes

Version	Section	Modifications
2020-09-25	All sections	Updates, addition of new sections Assay Procedure, History of Changes, Notes

 Manufacturer	 Date of manufacture	 Use by	LOT Batch code	REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	 Contains sufficient for <n> tests		